

Table of contents:

1. Logging onto the new Bruker 400 workstation	pg 1
2. How to change your password	pg 1
3. How to setup a basic new experiment	pg 2
4. Accessing your data	pg 4
4.1 Example	pg 5

1. Logging onto the new Bruker 400 workstation:

Computer name = **ucsf400.ucsf.edu**

Username = **contact the curator for this**


Password = **contact the curator for this**

2. How to change your password:

- 1) Open a new terminal from the main Linux window and type "**passwd**" at the prompt.
- 2) You will be asked to enter your old password, then your new password (twice).

3. How to Setup a Basic New Experiment:

- 1) Go to the **Start** tab.
- 2) Left-Click on the "**Create Dataset**" Icon.



This creates a "**New**" experiment window.

3) Fill in slots 1 through 8, then click **OK** when done.

NOTES:

- A) Slot #1 (Experiment Name): Please do not include spaces.
- B) Slot #2 (Experiment Number): Iterate numbers to store multiple experiments under same folder. See the example in Section 4.1 below. Enter numbers only. No characters or symbols allowed.
- C) Slot #7 (Data Directory) should be set to: **/home/data/OWNER/operator**
- D) Slot #8 (Title/annotations): Annotate your expt. here – this is good for bookkeeping purposes.

The screenshot shows the 'New...' dialog box in the software. The fields and their values are as follows:

- NAME:** MY_SAMPLE (Slot 1)
- EXPNO:** 1 (Slot 2)
- PROCNO:** 1 (Slot 3)
- Experiment:** PROTON (Slot 4)
- Solvent:** H2O+D2O (Slot 5)
- DIR:** /home/data/jtgrp/Jordan (Slot 7)
- TITLE:** (Empty field, Slot 8)

Other options include 'Use current parameters' (unchecked), 'Execute 'getprosol'' (checked), and 'Keep parameters' (unchecked). The 'Receivers' field is set to 1. Buttons for 'OK', 'Cancel', 'More Info...', and 'Help' are at the bottom.

4. Accessing your data:

- You can transfer your data to your local network using the secure ftp program of your choice (e.g. WinSCP for PCs or Cyberduck for Macs). Computer name = **ucsf400.ucsf.edu**

NOTES:

- A) Access must be done through a secured protocol such as **sftp** or **scp** (port 22).
 B) Use your linux account username and password.

- The root directories are located here:

/home/users/OWNER

- The main data directories are located here:

/home/data/OWNER

- A suggested account structure in Bruker Topspin 3.5 is as follows:

Owner	(linux account)
→ Operator	(individual user)
→ Experiment name	(name of your experiment queue)
→ Experiment #	(experiment #)
→ pdata/Process #	(processed experiment #)

The **Operator**, **Experiment name**, **Experiment #**, and **Process#** are automatically created when you launch a **Create Dataset** job.

Owner = user name, linux account.

Operator = individual user or project name.

4.1 An example:

Owner = **my_lab**
Operator = **john-doe**
Experiment name = **my_sample**

Experiment # **1** (example = 1-D PROTON)

Experiment # **2** (example = 1-D CARBON)

→ The resulting data will be located in:

/home/data/**my_lab/john-doe/my_sample/1**/fid (PROTON)
/home/data/**my_lab/john-doe/my_sample/2**/fid (CARBON)

The number after **my_sample** corresponds to the Experiment #.

2-D experiments will be named "**ser**".

If you have processed your data using TopSpin on the workstation:

/home/data/**my_lab/john-doe/my_sample/1**/pdata (PROTON)
/home/data/**my_lab/john-doe/my_sample/2**/pdata (CARBON)

Under the **pdata** directory:

.../pdata/1/1r (Your processed data)
.../pdata/1/title (Text file containing your annotations, if any)

The number after **pdata** corresponds to the Process #.

Processed 2D experiments will be named "**2rr**".