


# Gateway Reading Frame A

([http://tools.invitrogen.com/content/sfs/vectors/rfa\\_seq.txt](http://tools.invitrogen.com/content/sfs/vectors/rfa_seq.txt))



**Plasmid name:** HSP68-lacZ

**Plasmid size:** 7516 bp

**Constructed by:** Kothary et al., Development 105:707

**Construction date:** 1989

**Comments/References:** Unique sites in polylinkers: XhoI, Sall, HindIII, PstI----  
NotI----Sall cuts twice, but can be used to purify insert; partial list of non-cutting  
enzymes: BglII, RsrII, SexA1, NsiI