

Enhanced Middle Molecule Clearance by a Biomimetic Dialyzer Membrane

J.J. Groszek¹, J. Cheng¹, C. Blaha², R. Kant², J. Park², B. Chui³, K. Goldman⁴, S. Roy², W.H. Fissell^{1,5}

¹*Division of Nephrology and Hypertension, Vanderbilt University*

²*Department of Bioengineering and Therapeutic Sciences, UCSF*

³*Ben Chui Consulting*

⁴*H-Cubed, Inc.*

⁵*Department of Biomedical Engineering, Vanderbilt University*

Background:

Although polymer dialyzers attain very high small solute clearance rates, polydisperse pore sizes limit middle molecule clearance by conventional polymer dialyzers. We hypothesized that a biomimetic membrane with uniform slit pores would enhance middle molecule clearance.

Methods:

Thin film silicon membrane samples (surface area $1.4 \times 10^{-5} \text{ m}^2$) with monodisperse slit-shaped 5-8 nanometer pores ("SNM-HD") were manufactured as previously described, surface-modified with polyethylene glycol and mounted in a custom designed cartridge. Membrane pore size was estimated from hydraulic permeability measurements. Human blood was spiked with PABA and fluorescent-labelled Ficoll, a polydisperse globular polysaccharide. Blood flow was set at 100 ml/min and dialysate flow varied between 70 and 140 ul/min. Blood and dialysate concentrations of PABA and Ficoll were measured by size-exclusion chromatography. Clearance as a function of dialysate flow rate was extrapolated to estimate KoA as a function of molecular weight.

Results:

B2M clearance in high-flux dialyzers is typically 5-8% of the value of urea clearance. PABA clearance was 139 ml/min/m². 15.7Angstrom Ficoll (same radius as B2M) clearance in the SNM dialyzers was 16.9 ml/min/m², 3-5 fold higher than reported B2M clearance in polymer dialyzers.

Conclusions:

A biomimetic membrane of uniform slit pores potentially offers much higher middle molecule clearance than conventional polymer dialyzers.